

Competence creates Confidence.





• Model no. H3026

COOLING CHAMBER FOR CONDITIONING SPECIMENS

ISO 3126

The walk-in cooling chamber is used to condition specimens to the preset temperature.



Standard features

- | | |
|--|---------------------------------------|
| ● Access via one or two swing doors (depending on version) | ● Interior lighting with light switch |
| ● High temperature accuracy | ● Unit with compressor |
| ● Prefabricated cooling unit, air-cooled | ● Polyurethane foam wall insulation |
| ● Stainless steel base | |

Version COOLING CHAMBER FOR CONDITIONING SPECIMENS

		H3026-0021	H3026-0022	H3026-0023
Temperature range	°C	0 to -20		
Thermal insulation of the cell elements		Rigid polyurethane foam foamed with cyclopentane		
Wall thickness	mm	100		
Foam density thermal insulation	kg/m ³	40		
U-value	W/m ² K	0.20 / according to DIN EN 13165		
Frigen refrigerant		452 A		
Colour		White		
Base design		KR11 - Stainless steel base (1.4301)		
Slip resistance class	R	11		
Permissible wheel load	N	1000	4000	4000
Load area	cm ²	> 4		
Wheel material		Rubber		
Permissible surface load	N/m ²	30000	50000	50000
Ventilation of the underneath of the cooling chamber	m ²	1.44	5.67	5.67
Internal dimensions (length x width x height)	mm	1000 x 1000 x 1980	2500 x 1900 x 2280	2700 x 2500 x 2730
External dimensions (length x width x height)	mm	1200 x 1200 x 2180	2700 x 2100 x 2480	2900 x 2700 x 2930
Voltage data		400 V, 50 Hz Special voltages on request		